Instructions for use / Notice

Actifish

Please read carefully these instructions before using this product. They contain important information for the processing of fish semen.

DESCRIPTION OF THE PRODUCT

Actifish is a concentrated isotonic solution that activates spermatozoa and improves milt fertility.

Composition

Contains ultra-pure water, activation agent, salts, gentamycin

Antibiotic

Yes, 0.05 g/l in the concentrated solution

Technical characteristics

Concentrated (10X) solution for activation of fish spermatozoa to be prepared as indicated. Manufactured by the ISO 9001:2008- certified laboratories of IMV Technologies.

Expiration

24 months from date of manufacture (manufacturing batch No. and date indicated on the label).

Sterilization

No

Color

Transparent

USES OF ACTIFISH:

1) Activation

In some fish, spermatozoa remains static in the seminal liquid. It can represent approximately 80% of the milt volume. Milt becomes mobile in the presence of Actifish. Actifish has been developed by the **French Institute of Agronomic Research (INRA)** in collaboration with the **French Syndicate of Breeders in Aquaculture (SYSAAF)** and IMV Technologies to provide hatcheries standard and optimized procedures for fish fertility (1997).

Based on the studies of Billard et al. (1977, INRA, France), Actifish significantly improves eyed eggs rate in rainbow trout (Haffray et al., 2008) and Atlantic salmon (unpublished data).

2) Improvement

Milt with low concentration or low motility gives lower fertility rates in salmonids hatchery. This low quality of the milt can come from:

- The fish
 - Due to gonad maturity collected too early or too late
 - Due to sexual inversion (sex-reversed males/neomales)
 - Due to low fertility performance of a genetic strain
 - Due to sub-optimal diet during maturation
- The milt
 - Conserved for several hours/days
 - Cryopreserved

3) Standardization

The quality of the activation solution will impact the fertilization rate and the eyed-eggs rate. It is therefore important to use a standard and controlled solution, independent of water quality.

Actifish

INSTRUCTIONS FOR USE

Preparation of Actifish

- Before use, the Actifish bottle has to be diluted in distilled water, demineralized, de-ionized water or equivalent quality.
- **Do not use mineral water, circuit water or tap water.**
- **The water does not have to be sterile.**
- **Do not use the concentrated solution without diluting.**
- Dilution ratios

	Salmonids	Marine Fish	Cyprinidae	Siluridae
Examples	<i>Salmo salar</i> <i>Oncorhynchus kisutch</i> <i>Salmo trutta</i> <i>Oncorhynchus mykiss</i> <i>Oncorhynchus tshawytscha</i> <i>Salvelinus alpinus</i> <i>Salvelinus fontinalis</i> Other salmonids...	<i>Dicentrarchus labrax</i> <i>Sparus aurata</i> <i>Argyrosomus regius</i> <i>Lates calcarifer</i> Other marine fish...	<i>Carassius carassius</i> <i>Cyprinus carpio</i> <i>Labeo rohita</i> <i>Hypophthalmichthys molitrix</i> Other Cyprinidae...	<i>Silurus</i> Other...
Concentrated Actifish	1 vol	1 vol	1 vol	1 vol
Distilled Water	9 vol	3 vol	29 vol	29 vol
Total	10 vol	4 vol	30 vol	30 vol

- After preparation, the solution can be stored for 1 week at 4°C.

Ovocytes preparation

- After stripping of the ovocytes from the female, rinse them with Ovafish (Ref. 026950). See Ovafish instructions for detailed protocol.
- Place the ovocytes in a dry and clean bucket.

Milt preparation


- When the volume of milt is small, it is advised to first **dilute milt with Storfish** (Ref. 018500) to maximize the chances of spermatozoa to reach the ovocytes. See Storfish instructions for detailed protocol.
- If the milt is not stripped but collected from gonad dissection, it needs to **first be matured in Storfish** (Ref. 018500). See Storfish instructions for detailed protocol.

Fertilization

- Place the ovocytes in a dry and clean bucket.
- Add diluted milt to reach a ratio of 0.5-2 Million spermatozoa / ovocyte.
- Mix gently with the hand during 10 s to homogeneously distribute the spermatozoa in the entire bucket.

⚠ Use non-spermicide gloves to mix. Some plastic can release powder or toxic compounds dangerous for the cells. We recommend the Copolymer soft-touch glove (Ref. 024183).

- Add diluted Actifish 10 cm (about 10 fingers width) above the ovocyte level.

Instructions for use / Notice


Actifish

⚠ Avoid temperature shock: the Actifish must be approximately as the same temperature as the ovocytes ($\pm 5^{\circ}\text{C}$).

- Mix gently during 30 s.
- Wait for 2 minutes
- Remove the Actifish : **remove Actifish surplus**
- Disinfect the eggs if necessary and incubate

Sperm motility estimate

- Actifish is also adapted to standardize sperm quality estimates. The full protocol of fish milt quality control is available with IMV CASA system (Computer-Assisted-Sperm-Analysis, IVOS II, Ref.24911 and CEROS II, Ref.24905).
- Predilute the semen with Storfish in order to reach a concentration of approximately 200 Million spermatozoid / ml.
- Fill an Eppendorf tube with 1 ml of diluted Actifish. A few seconds before the analysis **for 100 μl** of prediluted semen in the Eppendorf cap

- Prepare the pipette to load 2.7 μl in Leja chamber (20 μm depth, Ref. 025107)
- Close the cap mix semen and diluted Actifish.
- Take the 2.7 μl of activated semen and rapidly load the chamber under the microscope and start acquisition after stabilization of the sample drift

STORAGE

- Do not use in case of visible deterioration of the packaging or the product, such as particles floating or color change.
- Store at temperatures between $+2^{\circ}\text{C}$ and $+8^{\circ}\text{C}$. Keep in a cool, dry, dark place. Avoid sunlight.
- Do not use after freezing.
- Do not use after the expiration date.
- Do not store close to foodstuff, drinks or animal feed. Do not ingest or inhale. In case of accidental projection on skin or eyes, rinse profusely with water. Consult a doctor in case of problem.

TRANSPORT

- Can be transported between 4°C and 35°C .

PRESENTATION

Description	Packaging	Reference
Actifish	Box of 10 x 1 bottle 1L of concentrated solution	018274

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