## Instructions for use / Notice

# FreezeFish

Please read carefully these instructions before using this product. They contain important information for the processing of **fish** milt.

### DESCRIPTION OF THE PRODUCT

**FreezeFish** has been developed specifically to improve milt cryopreservation of **freshwater fish**.

#### Composition

Contains ultra-pure water, salts, biological buffer, DMSO and lipid solution.

#### Antibiotics

No

#### Technical characteristics

Cryopreservation medium **to be prepared** as indicated and manufactured by the ISO 9001:2008- certified laboratories of IMV Technologies.

#### Expiration

12 months from date of manufacture (manufacturing batch No. and date indicated on the label).

#### Sterilization

Sterile by filtration

#### Color

Clear

### INSTRUCTIONS FOR USE

- If the DMSO is solid, take the DMSO out of the fridge. **After one hour at room temperature**, DMSO is liquid.
- To prepare FreezeFish media, open the kit and :
  - Open the bottle of Freezesol (437.5 ml)
  - Pour into it the content of the bottle of Freezelip (12.5ml)
  - Pour into it the bottle of DMSO (50 ml)
  - Mix gently
- After preparation, the solution can be stored for 1 week at 4°C
- Pour 1 volume of milt in 3 volumes of FreezeFish medium.
  - ⚠ **Avoid temperature shock:** the mix must be approximately as the same temperature as the milt ( $\pm 5^{\circ}\text{C}$ )
- Keep the preparation at 4°C (on ice or in the fridge) prior to freezing, or freeze it immediately. **Avoid keeping the preparation for more than 30 minutes** before freezing.
- Freeze with IMV protocol for freshwater fish using straws (Ref: 005569, 005702, 005703), automated filling and sealing machine MRS4 (Ref: 019995) and Digitcool (Ref: 007262) programmable freezer.
  - ⚠ If the milt is not stripped but collected from gonad dissection, it needs to **first be matured in Storfish** (Ref. 018500). See Storfish instructions for detailed protocol.

## Instructions for use / Notice

### FreezeFish

#### STORAGE

- Do not use in case of visible deterioration of the packaging or the product, such as particles floating.
- Store at temperatures between +2°C and +8°C. Keep in a cool, dry, dark place. Avoid sunlight.
- Do not use after freezing
- Do not use after the expiration date.
- Do not store close to foodstuff, drinks or animal feed. Do not ingest or inhale. In case of accidental projection on skin or eyes, rinse profusely with water. Consult a doctor in case of problem.

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