Instructions for use / Notice

Ovafish

Please read carefully these instructions before using this product. They contain important information for the processing of fish semen.

DESCRIPTION OF THE PRODUCT

Ovafish is a concentrated solution that washes salmonids ovocytes before fertilization

Composition

Contains ultra-pure water, salts and a biological buffer

Antibiotics

No

Technical characteristics

Concentrated (10X) solution for ovocytes washing to be prepared as indicated. Manufactured by the ISO 9001:2008-certified laboratories of IMV Technologies.

Expiration

24 months from date of manufacture (manufacturing batch No. and date indicated on the label).

Sterilization

Sterile by filtration

Color

Transparent

USES OF OVAFISH

Ovafish has been developed by the French institute of Agronomic research in collaboration with the French Syndicate of Breeders in Aquaculture and IMV Technologies to provide French breeders standard and optimized procedures for fish fertility (1997). Ovafish is a rinsing solution to wash ovocytes (Bonnet et al. (2003), Rime et al. (2004). It allows to remove faeces, blood, mucus and coelomic fluid and provides a significant increase in eyed eggs rate in rainbow trout (Haffray et al., 2008).

Ovafish is also used to conserve ovocytes in a clean and standard environment. Decreasing bacterial and viral loads in the ovocytes, Ovafish also limits the risk of vertically-transmitted disease (unpublished data). Use of Ovafish is recommended for salmonids, catfish, sturgeon and tilapia.

INSTRUCTIONS FOR USE

Preparation of Ovafish

- Before use, the Ovafish bottle has to be diluted in distilled water, demineralized, de-ionized water or equivalent quality.
- **Do not use mineral water, circuit water or tap water.**
- **Do not use the concentrated solution without diluting.**
- **The water does not have to be sterile.**
- The dilution ratio is 1 volume of concentrated Ovafish diluted in 9 volumes of distilled water in order to obtain 10 volumes of solution ready to use.
- After preparation, the solution can be stored for 1 week at 4°C

Ovocytes rinsing

- Place some diluted Ovafish in a dry and clean bucket.
- Strip the ovocytes into the bucket


⚠ Avoid temperature shock: Ovafish must be approximately **as the same temperature as the ovocytes** ($\pm 5^{\circ}\text{C}$).

Instructions for use / Notice

Ovafish

- Add diluted Ovafish up to 2 cm (about 2 fingers width) above the ovocytes level.
- Mix gently during 30 s to clean the ovocytes.
- ⚠ Use non-spermicide gloves to mix. Some plastic can release powder or toxic compounds dangerous for the cells. We recommend the copolymer soft-touch glove (Ref. **024183**).
- Remove the ovocytes with a strainer

Storage of ovocytes

- Place some diluted Ovafish in a dry and clean bucket.
- After rinsing the ovocytes as described previously, pour them in the bucket.
- Add diluted Ovafish up to 2 cm (about 2 fingers width) above the ovocytes level.
- **Store the bucket at 4°C and use it within one week**

Fertilization

- Remove the Ovafish with a strainer before fertilization
- Place the ovocytes in a dry and clean bucket.
- Add diluted milt and Actifish as described in Instructions of Storfish (Ref. 018500) and Actifish (Ref. 018274)

STORAGE

- Do not use in case of visible deterioration of the packaging or the product, such as particles floating or color change.
- Store at temperatures between +2°C and +8°C. Keep in a cool, dry, dark place. Avoid sunlight.
- Do not use after freezing.
- Do not use after the expiration date.
- Do not store close to foodstuff, drinks or animal feed. Do not ingest or inhale. In case of accidental projection on skin or eyes, rinse profusely with water. Consult a doctor in case of problem.

TRANSPORT

- Can be transported between 4 °C and 35 °C.

PRESENTATION

Description	Packaging	Reference
Ovafish	Bottle of 1L	026950

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